

Recombinant Human AAAS

Catalog No: #GP10972



Package Size: #GP10972-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human AAAS
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human achalasia, adrenocortical insufficiency, alacrimia
Target Name	achalasia, adrenocortical insufficiency, alacrimia
Other Names	AAA; AAASb; GL003; ALADIN; ADRACALA; ADRACALIN
Accession No.	Swissprot:Q9NRG9Gene Accession:BC000659
Uniprot	Q9NRG9
GeneID	8086;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a member of the WD-repeat family of regulatory proteins and may be involved in normal development of the peripheral and central nervous system. The encoded protein is part of the nuclear pore complex and is anchored there by NDC1. Defects in this gene are a cause of achalasia-addisonianism-alacrima syndrome (AAAS), also called triple-A syndrome or Allgrove syndrome. Two transcript variants encoding different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only