

Recombinant Human MLLT10

Catalog No: #GP10974



Package Size: #GP10974-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human MLLT10
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	myeloid/lymphoid or mixed-lineage leukemia; translocated to, 10
Other Names	AF10
Accession No.	Swissprot:P55197Gene Accession:BC080577
Uniprot	P55197
GeneID	8028;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a transcription factor and has been identified as a partner gene involved in several chromosomal rearrangements resulting in various leukemias. Multiple transcript variants encoding different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only