

Recombinant Human MLLT3

Catalog No: #GP10975



Package Size: #GP10975-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human MLLT3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 322-568 amino acids of human MLLT3
Target Name	MLLT3, super elongation complex subunit
Other Names	AF9; YEATS3
Accession No.	Swissprot:P42568Gene Accession:BC036089
Uniprot	P42568
GeneID	4300;
Storage	-20~-80°C, pH 7.6 PBS

Background

Component of the super elongation complex (SEC), a complex required to increase the catalytic rate of RNA polymerase II transcription by suppressing transient pausing by the polymerase at multiple sites along the DNA.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only