

Recombinant Human AGBL2

Catalog No: #GP10978



Package Size: #GP10978-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human AGBL2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 618-902 amino acids of human AGBL2
Target Name	ATP/GTP binding protein-like 2
Other Names	CCP2
Accession No.	Swissprot:Q5U5Z8Gene Accession:BC028200
Uniprot	Q5U5Z8
GeneID	79841;
Storage	-20~-80°C, pH 7.6 PBS

Background

AGBL2 is a 902 amino acid cytoplasmic protein belonging to the peptidase M14 family. CCP2 is considered a metalloprotease that may play a role in the processing of tubulin. CCP2 binds one zinc ion per subunit as a cofactor and exists as three alternatively spliced isoforms.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only