

Recombinant Human AGL5

Catalog No: #GP10979



Package Size: #GP10979-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human AGL5
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 300 amino acids of human ATP/GTP binding protein-like 5
Target Name	ATP/GTP binding protein-like 5
Accession No.	Swissprot:Q8NDL9Gene Accession:BC007415
Uniprot	Q8NDL9
GeneID	60509;
Storage	-20~-80°C, pH 7.6 PBS

Background

Metalloprotease that mediates tubulin deglutamylation. Specifically catalyzes the deglutamylation of the branching point glutamate side chains generated by post-translational glutamylation in proteins such as tubulins. In contrast, it is not able to act as a long-chain deglutamylase that shortens long polyglutamate chains, a process catalyzed by AGTPBP1/ CCP1, AGL1/CCP4 and AGL4/CCP6.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only