

Recombinant Human AGFG1

Catalog No: #GP10981



Package Size: #GP10981-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human AGFG1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250amino acids of human ArfGAP with FG repeats 1
Target Name	ArfGAP with FG repeats 1
Other Names	HRB; RAB; RIP
Accession No.	Swissprot:P52594Gene Accession:BC096273
Uniprot	P52594
GeneID	3267;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is related to nucleoporins, a class of proteins that mediate nucleocytoplasmic transport. The encoded protein binds the activation domain of the human immunodeficiency virus Rev protein when Rev is assembled onto its RNA target, and is required for the nuclear export of Rev-directed RNAs. Multiple transcript variants encoding different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only