

Recombinant Human AGPAT4

Catalog No: #GP10983



Package Size: #GP10983-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human AGPAT4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 146-306 amino acids of human 1-acylglycerol-3-phosphate O-acyltransferase 4
Target Name	1-acylglycerol-3-phosphate O-acyltransferase 4
Other Names	1-AGPAT4; dJ473J16.2; LPAAT-delta; RP3-473J16.2
Accession No.	Swissprot:Q9NRZ5Gene Accession:BC020209
Uniprot	Q9NRZ5
GeneID	56895;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the 1-acylglycerol-3-phosphate O-acyltransferase family. This integral membrane protein converts lysophosphatidic acid to phosphatidic acid, the second step in de novo phospholipid biosynthesis.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only