

Recombinant Human AGTRAP

Catalog No: #GP10986



Package Size: #GP10986-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human AGTRAP
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 77-159 amino acids of human angiotensin II receptor-associated protein
Target Name	angiotensin II receptor-associated protein
Other Names	ATRAP
Accession No.	Swissprot:Q6RW13Gene Accession:BC017328
Uniprot	Q6RW13
GenID	57085;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a transmembrane protein localized to the plasma membrane and perinuclear vesicular structures. The gene product interacts with the angiotensin II type I receptor and negatively regulates angiotensin II signaling. Alternative splicing of this gene generates multiple transcript variants encoding different isoforms.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only