

NUDT18 Antibody

Catalog No: #42946

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NUDT18 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total NUDT18 protein.
Immunogen Description	Fusion protein of human NUDT18
Target Name	NUDT18
Other Names	MTH3
Accession No.	Swiss-Prot#: Q6ZVK8Gene ID: 79873
Uniprot	Q6ZVK8
GeneID	79873;
Concentration	1.9mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Lung cancer tissue using #42946 at dilution 1/100,

Background

The protein encoded by this gene is a member of the Nudix hydrolase family. Nudix hydrolases eliminate potentially toxic nucleotide metabolites from the cell and regulate the concentrations and availability of many different nucleotide substrates, cofactors, and signaling molecules. This protein contains a Nudix hydrolase domain and hydrolyzes oxidized forms of guanosine and deoxyguanosine diphosphates.?

Note: This product is for in vitro research use only