

Recombinant Human AKR1B10

Catalog No: #GP10992



Package Size: #GP10992-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human AKR1B10
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	aldo-keto reductase family 1, member B10 (aldose reductase)
Other Names	HIS; HSI; ARL1; ARL-1; ALDRLn; AKR1B11; AKR1B12
Accession No.	Swissprot:O60218Gene Accession:BC008837
Uniprot	O60218
GeneID	57016;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the aldo/keto reductase superfamily, which consists of more than 40 known enzymes and proteins. This member can efficiently reduce aliphatic and aromatic aldehydes, and it is less active on hexoses. It is highly expressed in adrenal gland, small intestine, and colon, and may play an important role in liver carcinogenesis.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only