

## Recombinant Human ALG9

Catalog No: #GP11006



Package Size: #GP11006-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ALG9
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 427-611 amino acids of human ALG9
Target Name	ALG9, alpha-1,2-mannosyltransferase
Other Names	CDG1L; DIBD1; LOH11CR1J
Accession No.	Swissprot:Q9H6U8Gene Accession:BC009255
Uniprot	Q9H6U8
GeneID	79796;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes an alpha-1,2-mannosyltransferase enzyme that functions in lipid-linked oligosaccharide assembly. Mutations in this gene result in congenital disorder of glycosylation type II. Multiple transcript variants encoding different isoforms have been found for this gene.?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only