

Recombinant Human ORM1

Catalog No: #GP11010



Package Size: #GP11010-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ORM1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 19-201 amino acids of human orosomuroid 1
Target Name	orosomuroid 1
Other Names	ORM; AGP1; AGP-A
Accession No.	Swissprot:P02763Gene Accession:BC026238
Uniprot	P02763
GeneID	5004;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a key acute phase plasma protein. Because of its increase due to acute inflammation, this protein is classified as an acute-phase reactant. The specific function of this protein has not yet been determined; however, it may be involved in aspects of immunosuppression.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only