

Recombinant Human CRYAA

Catalog No: #GP11015



Package Size: #GP11015-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CRYAA
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-173 amino acids of human CRYAA
Target Name	alpha-crystallin A chain
Accession No.	Swissprot:P02489Gene Accession:BC069528
Uniprot	P02489
GeneID	102724652;1409;
Storage	-20~-80°C, pH 7.6 PBS

Background

Contributes to the transparency and refractive index of the lens. Has chaperone-like activity, preventing aggregation of various proteins under a wide range of stress conditions.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only