

# Recombinant Human ALPK1

Catalog No: #GP11019



Package Size: #GP11019-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ALPK1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein of human ALPK1
Target Name	alpha kinase 1
Other Names	LAK; 8430410J10Rik
Accession No.	Swissprot:Q96QP1 Gene Accession:BC060780
Uniprot	Q96QP1
GeneID	80216;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes an alpha kinase. Mice which were homozygous for disrupted copies of this gene exhibited coordination defects (PMID: 21208416). Multiple transcript variants encoding different isoforms have been found for this gene.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only