

Recombinant Human NIF3L1

Catalog No: #GP11021



Package Size: #GP11021-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human NIF3L1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 247-377 amino acids of human NIF3 NGG1 interacting factor 3-like 1 (<i>S. cerevisiae</i>)
Target Name	NIF3 NGG1 interacting factor 3-like 1 (<i>S. cerevisiae</i>)
Other Names	CALS-7; MDS015; ALS2CR1
Accession No.	Swissprot:Q9GZT8Gene Accession:BC007654
Uniprot	Q9GZT8
GeneID	60491;
Storage	-20~-80°C, pH 7.6 PBS

Background

NIF3-like protein 1 is a protein that in humans is encoded by the NIF3L1 gene. NIF3L1 has been shown to interact with COPS2. Expression analysis on cell lines revealed strong overexpression in a spermatogonia-derived cell line and a teratocarcinoma cell line. A human NIF3L1/GFP fusion protein showed localization of the protein in the cytoplasm.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only