

Recombinant Human ALS2CR11

Catalog No: #GP11024



Package Size: #GP11024-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ALS2CR11
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 103 amino acids of human amyotrophic lateral sclerosis 2 (juvenile) chromosome region, candidate 11
Target Name	amyotrophic lateral sclerosis 2 (juvenile) chromosome region, candidate 11
Accession No.	Swissprot:Q53TS8Gene Accession:BC030659
Uniprot	Q53TS8
GeneID	151254;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene has highly expression in digestive system component and liver, biliary system. It may have four isoforms that produced by alternative splicing. Transcription profiling of human mesenchymal stem cells transformed with lentiviral vectors of TERT, E6, SV40 and RAS.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only