

Recombinant Human AES

Catalog No: #GP11029



Package Size: #GP11029-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human AES
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 3-131 amino acids of human amino-terminal enhancer of split
Target Name	amino-terminal enhancer of split
Other Names	GRG; ESP1; GRG5; TLE5; AES-1; AES-2
Accession No.	Swissprot:Q08117Gene Accession:BC113735
Uniprot	Q08117
GenID	166;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is similar in sequence to the amino terminus of *Drosophila* enhancer of split groucho, a protein involved in neurogenesis during embryonic development. The encoded protein, which belongs to the groucho/TLE family of proteins, can function as a homooligomer or as a heterooligomer with other family members to dominantly repress the expression of other family member genes. Three transcript variants encoding different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only