

Recombinant Human AMOTL2

Catalog No: #GP11035



Package Size: #GP11035-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human AMOTL2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 200 amino acids of human angiotensin like 2
Target Name	angiotensin like 2
Other Names	LCCP
Accession No.	Swissprot:Q9Y2J4Gene Accession:BC011454
Uniprot	Q9Y2J4
GeneID	51421;
Storage	-20~-80°C, pH 7.6 PBS

Background

Angiotensin is a protein that binds angiotensin, a circulating inhibitor of the formation of new blood vessels (angiogenesis). Angiotensin mediates angiotensin inhibition of endothelial cell migration and tube formation in vitro. The protein encoded by this gene is related to angiotensin and is a member of the motins protein family.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only