

# Recombinant Human ANGPTL7

Catalog No: #GP11036



Package Size: #GP11036-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ANGPTL7
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 122-343 amino acids of human angiopoietin-like 7
Target Name	angiopoietin-like 7
Other Names	AngX; CDT6; dJ647M16.1; RP4-647M16.2
Accession No.	Swissprot:O43827Gene Accession:BC001881
Uniprot	O43827
GeneID	10218;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Angiopoietin-related protein 7 is a protein that in humans is encoded by the ANGPTL7 gene, a member of the ANGPTL family. Highly and specifically expressed in the cornea where is confined to the stromal layer.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only