

# Recombinant Human ANKRD54

Catalog No: #GP11042



Package Size: #GP11042-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ANKRD54
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 161 amino acids of human ankyrin repeat domain 54
Target Name	ankyrin repeat domain 54
Other Names	LIAR
Accession No.	Swissprot:Q6NXT1 Gene Accession:BC014641
Uniprot	Q6NXT1
GeneID	129138;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Plays an important role in regulating intracellular signaling events associated with erythroid terminal differentiation. Interacts (via ankyrin repeat region) with LYN (via SH3-domain) in an activation-independent status of LYN. By similarity. Forms a multiprotein complex with LYN and HCLS1. By similarity. Interacts with TSN2, VAV1, DBNL AND LASP1.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only