

Recombinant Human ANXA10

Catalog No: #GP11043



Package Size: #GP11043-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ANXA10
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	annexin A10
Other Names	ANX14
Accession No.	Swissprot:Q9UJ72?Gene Accession:BC007320
Uniprot	Q9UJ72
GeneID	11199;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. The function of this gene has not yet been determined.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only