

Recombinant Human AP2S1

Catalog No: #GP11047



Package Size: #GP11047-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human AP2S1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	adaptor-related protein complex 2, sigma 1 subunit
Other Names	AP17; HHC3; CLAPS2; AP17-DELTA
Accession No.	Swissprot:P53680Gene Accession:BC006337
Uniprot	P53680
GeneID	1175;
Storage	-20~-80°C, pH 7.6 PBS

Background

One of two major clathrin-associated adaptor complexes, AP-2, is a heterotetramer which is associated with the plasma membrane. This complex is composed of two large chains, a medium chain, and a small chain. This gene encodes the small chain of this complex. Alternative splicing has been observed in this gene and results in two known transcripts.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only