

# Recombinant Human ATG12

Catalog No: #GP11056



Package Size: #GP11056-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ATG12
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	autophagy related 12
Other Names	APG12; FBR93; APG12L; HAPG12
Accession No.	Swissprot:O94817Gene Accession:BC011033
Uniprot	O94817
GeneID	9140;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Autophagy is a process of bulk protein degradation in which cytoplasmic components, including organelles, are enclosed in double-membrane structures called autophagosomes and delivered to lysosomes or vacuoles for degradation. ATG12 is the human homolog of a yeast protein involved in autophagy.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only