

Recombinant Human APOC1

Catalog No: #GP11059



Package Size: #GP11059-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human APOC1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	apolipoprotein C-I
Accession No.	Swissprot:P02654Gene Accession:BC009698
Uniprot	P02654
GeneID	341;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a member of the apolipoprotein C1 family. This gene is expressed primarily in the liver, and it is activated when monocytes differentiate into macrophages. A pseudogene of this gene is located 4 kb downstream in the same orientation, on the same chromosome. This gene is mapped to chromosome 19, where it resides within a apolipoprotein gene cluster. Alternatively spliced transcript variants have been found for this gene, but the biological validity of some variants has not been determined.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only