

# Recombinant Human APOM

Catalog No: #GP11060



Package Size: #GP11060-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human APOM
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 23-188 amino acids of human apolipoprotein M
Target Name	apolipoprotein M
Other Names	G3a; NG20; apo-M; HSPC336
Accession No.	Swissprot:O95445Gene Accession:BC020683
Uniprot	O95445
GeneID	55937;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is an apolipoprotein and member of the lipocalin protein family. It is found associated with high density lipoproteins and to a lesser extent with low density lipoproteins and triglyceride-rich lipoproteins. The encoded protein is secreted through the plasma membrane but remains membrane-bound, where it is involved in lipid transport. Alternate splicing results in both coding and non-coding variants of this gene.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only