

Recombinant Human APOBEC4

Catalog No: #GP11061



Package Size: #GP11061-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human APOBEC4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 4 (putative)
Target Name	apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 4 (putative)
Other Names	C1orf169
Accession No.	Swissprot:Q8WW27Gene Accession:BC021711
Uniprot	Q8WW27
GenID	403314;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the AID/APOBEC family of polynucleotide (deoxy)cytidine deaminases, which convert cytidine to uridine. Other AID/APOBEC family members are involved in mRNA editing, somatic hypermutation and recombination of immunoglobulin genes, and innate immunity to retroviral infection.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only