

CNPY2 Antibody

Catalog No: #42954

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CNPY2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CNPY2 protein.
Immunogen Type	protein
Immunogen Description	Full length fusion protein of human CNPY2
Target Name	CNPY2
Other Names	MSAP; TMEM4; ZSIG9; HP10390
Accession No.	Swiss-Prot#: Q9Y2B0 Gene ID: 10330
Uniprot	Q9Y2B0
GeneID	10330;
Concentration	0.7mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Esophagus cancer tissue using #42954 at dilution 1/100.

Background

Positive regulator of neurite outgrowth by stabilizing myosin regulatory light chain (MRLC). It prevents MIR-mediated MRLC ubiquitination and its subsequent proteasomal degradation.

Note: This product is for in vitro research use only