

Recombinant Human ARHGAP15

Catalog No: #GP11076



Package Size: #GP11076-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ARHGAP15
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 280-470 amino acids of human Rho GTPase activating protein 15
Target Name	Rho GTPase activating protein 15
Other Names	BM046
Accession No.	Swissprot:Q53QZ3Gene Accession:BC016701
Uniprot	Q53QZ3
GenID	55843;
Storage	-20~-80°C, pH 7.6 PBS

Background

RHO GTPases regulate diverse biologic processes, and their activity is regulated by RHO GTPase-activating proteins (GAPs), such as ARHGAP15. GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state. Has activity toward RAC1. Overexpression results in an increase in actin stress fibers and cell contraction.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only