

## Recombinant Human ARHGEF2

Catalog No: #GP11078



Package Size: #GP11078-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ARHGEF2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 235-432 amino acids of human Rho/Rac guanine nucleotide exchange factor (GEF) 2
Target Name	Rho/Rac guanine nucleotide exchange factor (GEF) 2
Other Names	GEF; P40; GEFH1; LFP40; GEF-H1
Accession No.	Swissprot:Q92974Gene Accession:BC020567
Uniprot	Q92974
GenID	9181;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Rho GTPases play a fundamental role in numerous cellular processes that are initiated by extracellular stimuli that work through G protein coupled receptors. The encoded protein may form complex with G proteins and stimulate rho-dependent signals. Alternatively spliced transcript variants encoding different isoforms have been identified.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only