

Recombinant Human ARHGEF5

Catalog No: #GP11079



Package Size: #GP11079-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ARHGEF5
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human Rho guanine nucleotide exchange factor (GEF) 5
Target Name	Rho guanine nucleotide exchange factor (GEF) 5
Other Names	P60; TIM; GEF5; TIM1
Accession No.	Swissprot:Q12774Gene Accession:BC010046
Uniprot	Q12774
GenID	7984;
Storage	-20~-80°C, pH 7.6 PBS

Background

Rho GTPases play a fundamental role in numerous cellular processes initiated by extracellular stimuli that work through G protein coupled receptors. The encoded protein may form a complex with G proteins and stimulate Rho-dependent signals. This protein may be involved in the control of cytoskeletal organization.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only