

Recombinant Human ARHGEF1

Catalog No: #GP11081



Package Size: #GP11081-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ARHGEF1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 662-850 amino acids of human Rho guanine nucleotide exchange factor (GEF) 1
Target Name	Rho guanine nucleotide exchange factor (GEF) 1
Other Names	LSC; GEF1; LBCL2; SUB1.5; P115-RHOGEF
Accession No.	Swissprot:Q92888Gene Accession:BC034013
Uniprot	Q92888
GenID	9138;
Storage	-20~-80°C, pH 7.6 PBS

Background

Rho GTPases play a fundamental role in numerous cellular processes that are initiated by extracellular stimuli that work through G protein coupled receptors. The encoded protein may form complex with G proteins and stimulate Rho-dependent signals. Multiple alternatively spliced transcript variants have been found for this gene, but the full-length nature of some variants has not been defined.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only