

Recombinant Human DIRAS3

Catalog No: #GP11085



Package Size: #GP11085-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human DIRAS3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-226 amino acids of human DIRAS3
Target Name	DIRAS family GTPase 3
Other Names	ARHI; NOEY2
Accession No.	Swissprot:O95661Gene Accession:BC005362
Uniprot	O95661
GeneID	9077;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the ras superfamily. This gene is imprinted gene with monoallelic expression of the paternal allele which is associated with growth suppression. The encoded protein acts as a tumor suppressor whose function is abrogated in many ovarian and breast cancers. This protein may also play a role autophagy in certain cancer cells by regulating the autophagosome initiation complex.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only