

## Recombinant Human ARL4A

Catalog No: #GP11087



Package Size: #GP11087-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |  |
|-----------------------|--|
| Product Name          | Recombinant Human ARL4A  |
| Brief Description     | Recombinant Protein  |
| Immunogen Description | Fusion protein corresponding to a region derived from 2-197 amino acids of human ADP-ribosylation factor-like 4A |
| Target Name           | ADP-ribosylation factor-like 4A  |
| Other Names           | ARL4   |
| Accession No.         | Swissprot:P40617Gene Accession:BC001111  |
| Uniprot               | P40617   |
| GenID                 | 10124;   |
| Storage               | -20~-80°C, pH 7.6 PBS  |

## Background

ADP-ribosylation factor-like 4A is a member of the ADP-ribosylation factor family of GTP-binding proteins. ARL4A is similar to ARL4C and ARL4D and each has a nuclear localization signal and an unusually high guanine nucleotide exchange rate. ARL4A is located in both the nuclear and extranuclear cell compartments. Multiple transcript variants encoding the same protein have been found for this gene

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only