

# Recombinant Human ARL8A

Catalog No: #GP11088



Package Size: #GP11088-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ARL8A
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	ADP ribosylation factor like GTPase 8A
Other Names	GIE2; ARL10B
Accession No.	Swissprot:Q96BM9Gene Accession:BC015408
Uniprot	Q96BM9
GeneID	127829;
Storage	-20~-80°C, pH 7.6 PBS

## Background

May play a role in lysosomes motility. Alternatively, may play a role in chromosome segregation.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only