

# Recombinant Human ARMCX1

Catalog No: #GP11089



Package Size: #GP11089-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ARMCX1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human armadillo repeat containing, X-linked 1
Target Name	armadillo repeat containing, X-linked 1
Other Names	ALEX1
Accession No.	Swissprot:Q9P291Gene Accession:BC002691
Uniprot	Q9P291
GeneID	51309;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a member of the ALEX family of proteins and may play a role in tumor suppression. The encoded protein contains a potential N-terminal transmembrane domain and two Armadillo (arm) repeats. Other proteins containing the arm repeat are involved in development, maintenance of tissue integrity, and tumorigenesis. This gene is closely localized with other family members, including ALEX2 and ALEX3, on the X chromosome.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only