

Recombinant Human ACTR2

Catalog No: #GP11092



Package Size: #GP11092-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ACTR2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 300 amino acids of human ARP2 actin-related protein 2 homolog (yeast)
Target Name	ARP2 actin-related protein 2 homolog (yeast)
Other Names	ARP2
Accession No.	Swissprot:P61160Gene Accession:BC014546
Uniprot	P61160
GenID	10097;
Storage	-20~-80°C, pH 7.6 PBS

Background

The specific function of this gene has not yet been determined; however, the protein it encodes is known to be a major constituent of the ARP2/3 complex. This complex is located at the cell surface and is essential to cell shape and motility through lamellipodial actin assembly and protrusion. Two transcript variants encoding different isoforms have been found for this gene.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only