

# Recombinant Human GAB1

Catalog No: #GP11093



Package Size: #GP11093-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human GAB1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 149-449 amino acids of human GRB2-associated binding protein 1
Target Name	GRB2-associated binding protein 1
Accession No.	Swissprot:Q13480Gene Accession:BC064848
Uniprot	Q13480
GeneID	2549;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is a member of the IRS1-like multisubstrate docking protein family. It is an important mediator of branching tubulogenesis and plays a central role in cellular growth response, transformation and apoptosis. Two transcript variants encoding different isoforms have been found for this gene.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only