

Recombinant Human ARSA

Catalog No: #GP11094



Package Size: #GP11094-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | Recombinant Human ARSA |
| Brief Description | Recombinant Protein |
| Immunogen Description | Fusion protein corresponding to C terminal 250 amino acids of human arylsulfatase A |
| Target Name | arylsulfatase A |
| Other Names | MLD |
| Accession No. | Swissprot:P15289Gene Accession:BC014210 |
| Uniprot | P15289 |
| GeneID | 410; |
| Storage | -20~-80°C, pH 7.6 PBS |

Background

The protein encoded by this gene hydrolyzes cerebroside sulfate to cerebroside and sulfate. Defects in this gene lead to metachromatic leucodystrophy (MLD), a progressive demyelination disease which results in a variety of neurological symptoms and ultimately death. Alternatively spliced transcript variants have been described for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only