

Recombinant Human IFITM3

Catalog No: #GP11102

Package Size: #GP11102-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human IFITM3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to N terminal 107 amino acids of human interferon induced transmembrane protein 3
Target Name	interferon induced transmembrane protein 3
Other Names	1-8U; IP15; DSPA2b
Accession No.	Swissprot:Q01628Gene Accession:BC006794
Uniprot	Q01628
GenID	10410;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is an interferon-induced membrane protein that helps confer immunity to influenza A H1N1 virus, West Nile virus, and dengue virus. Two transcript variants, only one of them protein-coding, have been found for this gene. Another variant encoding an N-terminally truncated isoform has been reported, but the full-length nature of this variant has not been determined.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only