

# Recombinant Human PRDM14

Catalog No: #GP11104



Package Size: #GP11104-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human PRDM14
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 300 amino acids of human PR domain containing 14
Target Name	PR domain containing 14
Other Names	PFM11
Accession No.	Swissprot:Q9GZV8Gene Accession:BC052311
Uniprot	Q9GZV8
GeneID	63978;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a member of the PRDI-BF1 and RIZ homology domain containing (PRDM) family of transcriptional regulators. The encoded protein may possess histone methyltransferase activity and plays a critical role in cell pluripotency by suppressing the expression of differentiation marker genes. Expression of this gene may play a role in breast cancer

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only