

# Recombinant Human ARSK

Catalog No: #GP11107



Package Size: #GP11107-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ARSK
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 23-323 amino acids of human arylsulfatase family, member K
Target Name	arylsulfatase family, member K
Other Names	TSULF
Accession No.	Swissprot:Q6UWY0Gene Accession:BC130329
Uniprot	Q6UWY0
GeneID	153642;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Sulfatases (EC 3.1.5.6), such as ARSK, hydrolyze sulfate esters from sulfated steroids, carbohydrates, proteoglycans, and glycolipids. They are involved in hormone biosynthesis, modulation of cell signaling, and degradation of macromolecules.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only