

## Recombinant Human ART3

Catalog No: #GP11109



Package Size: #GP11109-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ART3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 27-362 amino acids of human ADP-ribosyltransferase 3
Target Name	ADP-ribosyltransferase 3
Other Names	ARTC3
Accession No.	Swissprot:Q13508Gene Accession:BC008397
Uniprot	Q13508
GeneID	419;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes an arginine-specific ADP-ribosyltransferase. The encoded protein catalyzes a reversible reaction which modifies proteins by the addition or removal of ADP-ribose to an arginine residue to regulate the function of the modified protein. An ADP-ribosyltransferase pseudogene is located on chromosome 11. Multiple transcript variants encoding different isoforms have been found for this gene.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only