

Recombinant Human SAE1

Catalog No: #GP11116



Package Size: #GP11116-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SAE1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	SUMO1 activating enzyme subunit 1
Other Names	AOS1; SUA1; UBLE1A; HSPC140
Accession No.	Swissprot:Q9UBE0Gene Accession:BC000344
Uniprot	Q9UBE0
GeneID	10055;
Storage	-20~-80°C, pH 7.6 PBS

Background

Posttranslational modification of proteins by the addition of the small protein SUMO, or sumoylation, regulates protein structure and intracellular localization. SAE1 and UBA2 form a heterodimer that functions as a SUMO-activating enzyme for the sumoylation of proteins.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only