

Recombinant Human BCL2L14

Catalog No: #GP11117



Package Size: #GP11117-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human BCL2L14
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to N terminal 253 amino acids of human BCL2-like 14 (apoptosis facilitator)
Target Name	BCL2-like 14 (apoptosis facilitator)
Other Names	BCLG
Accession No.	Swissprot:Q9BZR8Gene Accession:BC025778
Uniprot	Q9BZR8
GeneID	79370;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene belongs to the BCL2 protein family. BCL2 family members form hetero- or homodimers and act as anti- or pro-apoptotic regulators that are involved in a wide variety of cellular activities. Overexpression of this gene has been shown to induce apoptosis in cells. Three alternatively spliced transcript variants encoding two distinct isoforms have been reported for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only