

# Recombinant Human ALPI

Catalog No: #GP11121



Package Size: #GP11121-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ALPI
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 203-503 amino acids of human alkaline phosphatase, intestinal
Target Name	alkaline phosphatase, intestinal
Other Names	IAP
Accession No.	Swissprot:P09923Gene Accession:BC132678
Uniprot	P09923
GenID	248;
Storage	-20~-80°C, pH 7.6 PBS

## Background

There are at least four distinct but related alkaline phosphatases: intestinal, placental, placental-like, and liver/bone/kidney (tissue non-specific). The intestinal alkaline phosphatase gene encodes a digestive brush-border enzyme. This enzyme is upregulated during small intestinal epithelial cell differentiation.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only