

Recombinant Human RPL26L1

Catalog No: #GP11127



Package Size: #GP11127-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RPL26L1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	ribosomal protein L26-like 1
Other Names	RPL26P1
Accession No.	Swissprot:Q9UNX3Gene Accession:BC017360
Uniprot	Q9UNX3
GeneID	51121;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein that shares high sequence similarity with ribosomal protein L26. It is not currently known whether the encoded protein is a functional ribosomal protein or whether it has evolved a function that is independent of the ribosome. Transcript variants utilizing alternative polyA signals exist.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only