

Recombinant Human RPN1

Catalog No: #GP11133



Package Size: #GP11133-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RPN1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 458-607 amino acids of human ribophorin I
Target Name	ribophorin I
Other Names	OST1; RBPH1
Accession No.	Swissprot:P04843Gene Accession:BC010839
Uniprot	P04843
GeneID	6184;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a type I integral membrane protein found only in the rough endoplasmic reticulum. The encoded protein is part of an N-oligosaccharyl transferase complex that links high mannose oligosaccharides to asparagine residues found in the Asn-X-Ser/Thr consensus motif of nascent polypeptide chains. This protein forms part of the regulatory subunit of the 26S proteasome and may mediate binding of ubiquitin-like domains to this proteasome.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only