

Recombinant Human MTM1

Catalog No: #GP11148



Package Size: #GP11148-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human MTM1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 240-484 amino acids of human MTM1
Target Name	myotubularin 1
Other Names	CNM; MTMX; XLMTM
Accession No.	Swissprot:Q13496Gene Accession:BC030779
Uniprot	Q13496
GeneID	4534;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a dual-specificity phosphatase that acts on both phosphotyrosine and phosphoserine. It is required for muscle cell differentiation and mutations in this gene have been identified as being responsible for X-linked myotubular myopathy.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only