

# Recombinant Human PRPSAP1

Catalog No: #GP11155



Package Size: #GP11155-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human PRPSAP1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human phosphoribosyl pyrophosphate synthetase-associated protein 1
Target Name	phosphoribosyl pyrophosphate synthetase-associated protein 1
Other Names	PAP39
Accession No.	Swissprot:Q14558Gene Accession:BC009012
Uniprot	Q14558
GenelD	5635;
Storage	-20~-80°C, pH 7.6 PBS

## Background

PRPSAP1 (phosphoribosyl pyrophosphate synthetase-associated protein 1), also known as PAP39, is a 356 amino acid ubiquitous protein belonging to the ribose-phosphate pyrophosphokinase family. PRPSAP1 may play a negative regulatory role in 5-phosphoribose 1-diphosphate synthesis and is encoded by a gene mapping to human chromosome 17q25.1.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only