

Recombinant Human TWF2

Catalog No: #GP11162



Package Size: #GP11162-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TWF2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 177-349 amino acids of human twinfilin actin-binding protein 2
Target Name	twinfilin actin-binding protein 2
Other Names	A6r; A6RP; PTK9L; MSTP011
Accession No.	Swissprot:Q6IBS0Gene Accession:BC016452
Uniprot	Q6IBS0
GeneID	11344;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene was identified by its interaction with the catalytic domain of protein kinase C-zeta. The encoded protein contains an actin-binding site and an ATP-binding site. It is most closely related to twinfilin (PTK9), a conserved actin monomer-binding protein.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only